

Product Name :

Connecting Lanyard Fall Arrest Single Polyamide Lanyard NU4003

Product Code :

NLAB-SAFETYLABP20096



Description :

Connecting Lanyard Fall Arrest Single Polyamide Lanyard NU4003

Technical Specification :

Length: 2 metre from hook to hook.

Shock absorber with karabiner and other side scaffolding hooks.

Strength: 23 kN.

Rope dia: 12 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003