

**Product Name :**  
Connecting Lanyard Fall Arrest Single Polyamide Lanyard NU4004

**Product Code :**  
NLAB-SAFETYLABP20097



**Description :**

Connecting Lanyard Fall Arrest Single Polyamide Lanyard NU4004

**Technical Specification :**

Length: 2 metre from hook to hook.

Rope dia: 12 mm.

Strength: 23 kN.

Shock absorber with karabiner and other side scaffolding hooks.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003