

Product Name :
Connecting Lanyard Fall Arrest Single Webbing Lanyard NU7001

Product Code :
NLAB-SAFETYLABP200104



Description :

Connecting Lanyard Fall Arrest Single Webbing Lanyard NU7001

Technical Specification :

44 mm wide webbing energy absorbing single Lanyard.

23 kN strength.

Length: 2 metre from hook to hook.

Shock absorber with karabiner and other side scaffolding hooks.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003