

**Product Name :**  
Connecting Lanyard NJ 2004

**Product Code :**  
NLAB-SAFETYLABP200112



**Description :**

Connecting Lanyard NJ 2004

**Technical Specification :**

Work positioning braided rope lanyard for hands free work.

23 kN strength.

One side karabiner with aluminum rope adjuster and other side carabiner.

Works on 10.5 mm dia kermental rope.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003