

Product Name :
Connecting Lanyard NJ 3005

Product Code :
NLAB-SAFETYLABP200113



Description :

Connecting Lanyard NJ 3005

Technical Specification :

Works on 10.5 mm dia kermental rope.

One side karabiner with grip adjuster and other side carabiner.

Adjustable gradually and moves freely.

Can be accurately controlled with one hand.

New and compact grip device, used for hands free work.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003