

Product Name :
Connecting Lanyard Work Restraint Braided Lanyard NU 1004

Product Code :
NLAB-SAFETYLABP20082



Description :

Connecting Lanyard Work Restraint Braided Lanyard NU 1004

Technical Specification :

With one side karabiner and other side scaffolding hooks.

Strength: 23 kN.

Length: 2 metres from hook to hook.

Rope dia: 10.5 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003