

**Product Name :**  
Connecting Lanyard Work Restraint Braided Lanyard NU1001

**Product Code :**  
NLAB-SAFETYLABP20079



**Description :**

Connecting Lanyard Work Restraint Braided Lanyard NU1001

**Technical Specification :**

With one side eye loop and other side karabiner.

Length: 2 metres from hook to hook.

Rope dia: 10.5 mm.

Strength: 23 kN.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003