

**Product Name :**

Connecting Lanyard Work Restraint Elasticated Lanyard NU3003

**Product Code :**

NLAB-SAFETYLABP20091

**Description :**

Connecting Lanyard Work Restraint Elasticated Lanyard NU3003

**Technical Specification :**

With one side karabiner and other side hooks.

Length: Extendable from 1.2 metre to 2 metre.

Elasticated single Lanyard.

Extremely lightweight & compact lanyard.

35 mm wide 100% polyester webbing.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003