

Product Name :

Connecting Lanyard Work Restraint Elasticated Lanyard NU3004

Product Code :

NLAB-SAFETYLABP20092

**Description :**

Connecting Lanyard Work Restraint Elasticated Lanyard NU3004

Technical Specification :

With one side karabiner and other side scaffolding hook

Length: Extendable from 1.2 metre to 2 metre.

Elasticated single Lanyard.

Extremely lightweight & compact lanyard.

35 mm wide 100% polyester webbing.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003