

**Product Name :**  
Connecting Lanyard Work Restraint Elasticated Lanyard NU3005

**Product Code :**  
NLAB-SAFETYLABP20093



**Description :**

Connecting Lanyard Work Restraint Elasticated Lanyard NU3005

**Technical Specification :**

With one side karabiner and other side scaffolding hooks.

Elasticated single Lanyard.

Length: Extendable from 1.2 metre to 2 metre.

35 mm wide 100% polyester webbing.

Extremely lightweight & compact lanyard.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003