

Product Name :
Connecting Lanyard Work Restraint Webbing Lanyard NU2001

Product Code :
NLAB-SAFETYLABP20084



Description :

Connecting Lanyard Work Restraint Webbing Lanyard NU2001

Technical Specification :

With one side eye loop and other side karabiner.

Strength: 23 kN.

Length: 2 metres from hook to hook.

Webbing Width: 44 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003