

Product Name :
Consolidation App. (Three Gang)

Product Code :
NLAB-TESTINGLABM109042



Description :

Consolidation App. (Three Gang)

Technical Specification :

Rigidly manufactured from aluminous alloy casting to provide a high degree of accuracy with any frame distortions under load. The load bridge group is supported in high accuracy self-aligning seat balls. The beam provides three loading ratio: 9:1 10:1 11:1 and the beam assembly is fitted with an adjustable counterbalance weight. It consists of 3 consolidation test assemblies mounted on a single frame. Complete with 3 sets of Consolidation cell assembly for 60mm dia samples and three sets of weights each giving 10 kg/sq.cm. On each sample. And complete with rod holding the weights and coupling block holding the dial gauge or transducer. Supplied 'without': Consolidation cell, Weights, Dial Gauge (or transducer), holding bench which have to be ordered separately.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

