

Product Name :
Consolidation Test Apparatus

Product Code :
NLAB-CIVILLAB110001



Description :

Consolidation Test Apparatus

Technical Specification :

The apparatus has to be completed with the axial displacement gauges or transducers, Beam loading device(s) and Slotted weights. See accessories.

Overall dimensions (lxdxh): 1300x800x970 mm

Weight approx.: 145 kg

To perform simultaneously the anisotropic consolidation stage of three triaxial specimens

To reduce the testing time where only one triaxial compression machine is available

For triaxial cells up to 100 mm dia. specimens

Ergonomic design for better use of space in the laboratory

Suitable for other make of triaxial cells (with little adjustment)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003