

Product Name :
Constant Head Permeability Apparatus

Product Code :
NLAB-CIVILLAB110003



Description :

Constant Head Permeability Apparatus

Technical Specification :

The apparatus is used for the laboratory determination of permeability of soil using a constant or a variable head. This test is recommended for soils with co-efficient permeability in the range of 10^{-3} to 10^{-7} cm/sec. The maximum particle size of the soil which can be tested in the mould is 10mm. The soil permeability is a very important factor to study the behavior of soil in its natural condition with respect to water flow.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003