

Product Name :
Constant Pressure System for Rocks

Product Code :
NLAB-EDUSCHOOLEQ201001



Description :

Constant Pressure System for Rocks

Technical Specification :

The system consists of an oil pump, continuously driven by an electric motor during the entire period of operation to maintain the desired pressure.

The apparatus is designed to provide confining pressure upto 160 Bar to triaxial cells.

The unit provides continuous variable pressure upto 160 Bar, which can be increased or decreased, simply by turning a control knob.

Constant Pressure System for Rocks.

Constant Pressure System, 600 kg/cm².

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

