

**Product Name :**  
Continuous Stirred Tank Reactor India

**Product Code :**  
NLAB-VOCATIONALTVET29001



**Description :**

Continuous Stirred Tank Reactor

**Technical Specification :**

Feed Tank: Material Stainless Steel.

Capacity - 20 Ltrs.

Feed Circulation: By compressed air.

Flow Measurement: Rotameter 2Nos.

Pressure Regulator: 0-2 Kg/cm<sup>2</sup>

Reactor: Material Stainless Steel, Capacity 2 Ltrs (Approx).

Stirrer: Stainless Steel Impeller and shaft coupled with FHP Motor.

Pressure Gauge: Bourdon type 0-2 Kg/cm<sup>2</sup>

Stop Watch: Electronic

Control Panel comprises of Standard make the on-off switch, Mains Indicator etc.

Piping: Stainless Steel and PVC.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003