

Product Name : Controlled Drug Cabinet	Product Code : NLABMEDICAL-M3270009
Description :	
Controlled Drug Cabinet	
Technical Specification :	
<p>Medical drugs cabinets exceed the requirements of "The Misuse Of Drugs (Safe Custody) Regulations 1973" and fully comply with BS 2281.</p> <p>All cabinets are powder coated and oven baked for a tough, long lasting finish.</p> <p>All products are fitted with a 7 lever mortice dead lock.</p> <p>Minimum 16 gauge (approximately 1.6mm) mild steel used on all cabinets up to a maximum of 6mm on heavier duty cabinets.</p> <p>Floor or wall mounted options available.</p> <p>Optional alarm light and alarm with power supply available - factory fitted.</p> <p>Floor fixing kit available for selected cabinets.</p> <p>Multi-point locking featured on selected cabinets.</p> <p>CDC have partially adjustable shelves.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003