

Product Name : Controlled Drug Cabinet	Product Code : NLABMEDICAL-M3270009
--	---

--

Description :

Controlled Drug Cabinet

Technical Specification :

Medical drugs cabinets exceed the requirements of "The Misuse Of Drugs (Safe Custody) Regulations 1973" and fully comply with BS 2281.

All cabinets are powder coated and oven baked for a tough, long lasting finish.

All products are fitted with a 7 lever mortice dead lock.

Minimum 16 gauge (approximately 1.6mm) mild steel used on all cabinets up to a maximum of 6mm on heavier duty cabinets.

Floor or wall mounted options available.

Optional alarm light and alarm with power supply available - factory fitted.

Floor fixing kit available for selected cabinets.

Multi-point locking featured on selected cabinets.

CDC have partially adjustable shelves.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003