

Product Name : COP Tester or Oil Penetration Tester	Product Code : NLAB-TESTINGLABM41003
Description : COP Tester or Oil Penetration Tester	
Technical Specification : The apparatus is designed for determining the penetration time of oil or other liquid through paper and some other flat materials. A standard tilting type penetration testing apparatus consist of a chamber with Orifice 6.03cm diameter for holding the liquid medium (castor oil) and a clamping device for holding the paper. Digital temperature of 35 ± 10 C. Apparatus consists automatic digital timer for measuring the penetration time in seconds.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003