

Product Name :
Coplin Jar

Product Code :
NLAB-EDUSCHOOLEQ125001



Description :

Coplin Jar

Technical Specification :

Polylab Coplin Jar made of Polypropylene, can hold 10 slides of the size 1"x3" back to back.

The interior of the Coplin Jar is grooved in a way to hold slides in erect position.

Domed & shallow thread screw cap provides protection to the slides.

Packing

Coplin Jar

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003