

Product Name :
Cordless LED Binocular Stereo Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730108

Description :

Cordless LED Binocular Stereo Microscope

Technical Specification :

The LED illumination guarantees a clear and crisp, cool-to-touch white light, which lasts 60x longer than traditional tungsten bulbs.

This 45° inclined, binocular stereo microscope with diopter adjustment on one ocular tube, promises top viewing comfort.

The 2x objective is parfocaled and parcentered, and with a 10x WF and 20x WF eyepiece, offer magnifications of 20x and 40x.

The rechargeable battery enhances flexibility, adding a cordless feature, which guarantees easy transport and increased mobility.

The adjustable interpupillary distance ranges from 51mm to 75mm.

A dual illumination switch allows for three light setting options: top, bottom, or both.

Rugged cast aluminum with gray enamel finish.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003