

Product Name :
Corrosion Studies Kit

Product Code :
NLAB-TECHNICALAB37010



Description :

Corrosion Studies Kit

Technical Specification :

The corrosion studies kit allows simultaneous study of up to eight corrosion cells each containing three test specimens.

Sample mounting in a manner to minimise secondary effects.

Corrosion rates measured visually and by weighing.

Stirring by air or inert gas.

Comprising:

A service panel

Air pump

Eight test cells of glass construction, with specially machined lids enabling samples to be mounted

Test cell supports

Digital pH meter

Specimen pieces of steel, zinc, brass and copper

Platinum electrodes - Low voltage supply

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003