

Product Name :
Crack Measuring Microscope

Product Code :
NLAB -CONCRLAB120046



Description :

Crack Measuring Microscope

Technical Specification :

The apparatus operates by an adjustable lamp unit and the image is focused by turning a knob. The eyepiece scale can be turned through 360° to align with the direction of the crack or pitch under examination.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003