

Product Name :
CRC Code Decode Trainer

Product Code :
NLAB-ELECTRONICSAB160008



Description :

CRC Code Decode Trainer

Technical Specification :

Pulse Generator

One Switch to provide Low to High Pulse.

One Switch to provide High to Low Pulse.

Data Generator

8 bit Digital Data generator to generate any binary input word.

Selection of Data using 8-way DIP Switch.

Display

One LED indicator to indicate Power input.

On-board Circuits

CRC Encoder/ Transmitter Circuit for 16 bits CRC Polynomials

CRC Decoder/ Receiver Circuit for 16 bits CRC Polynomials.

Data Receiver Circuit for 8 Data bits.

Power Supply

Fixed DC power supply : + 5 V/500mA

Interconnections

All interconnections are made using 2mm banana Patch cords.

Bare board Tested Glass Epoxy PCB is used.

Set of 2mm Patch cords for interconnections.

User's Manual.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003