

Product Name :
Creep Testing Machine

Product Code :
NLAB-ENGINEERINGLB25004



Description :

Creep Testing Machine

Technical Specification :

Creep Testing Machine Lab Equipments Manufacturer, Suppliers & Exporters Specification of Creep Testing Machine

Table unit for the investigation of the creep behaviour of different materials

Flat specimens made of lead and PE

Specimen cross-section 2x5mm, measured length 25mm

Tensile stress range 5...25N/mm²

Position measurement gauge 0...10mm, resolution 0.01mm

Experiments above or below room temperature possible with conditioning box

Technical Data

Specimen

Material: lead, PE

Cross-section: 2x5mm

Measured length: 25mm

Tensile stress range: 5...25N/mm²

Position Measurement

Range: 0...10mm

Resolution 0.01mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003