

Product Name : CRO	Product Code : EC325
------------------------------	--------------------------------



Description :
CRO
Technical Specification :
<p>CRO</p> <p>4 Channel/2Ghz Digital Oscilloscope</p> <p>Digital real-time sampling: Accurately capture signals with at least 10X over-sampling on all channels, all the time with Tektronix proprietary sampling technology.</p> <p>Built-in analysis tools: Simplify analysis of your device with 16 automated measurements, FFT analysis, waveform math and cursors.</p> <p>Waveform limit testing: Eliminate mistakes with the pass/fail summary table that clearly shows how many waveforms are within the limits of your specified template.</p> <p>Help menu: Get the help you need “when and where you need it” with the built-in, context-sensitive help system.</p> <p>Dedicated front-panel controls: Spend less time learning and more time on your task with easy-to-use front-panel</p> <p>Affordable performance, compact design. 10x oversampling, USB connectivity, 16 automated measurements.</p> <p>Front-panel USB host port: Quickly store and transfer your waveforms and settings.</p>

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com