

Product Name : CRO	Product Code : EC325
------------------------------	--------------------------------



Description :

CRO

Technical Specification :

CRO
4 Channel/2Ghz Digital Oscilloscope
Digital real-time sampling: Accurately capture signals with at least 10X over-sampling on all channels, all the time with Tektronix proprietary sampling technology.
Built-in analysis tools: Simplify analysis of your device with 16 automated measurements, FFT analysis, waveform math and cursors.
Waveform limit testing: Eliminate mistakes with the pass/fail summary table that clearly shows how many waveforms are within the limits of your specified template.
Help menu: Get the help you need – when and where you need it – with the built-in, context-sensitive help system.
Dedicated front-panel controls: Spend less time learning and more time on your task with easy-to-use front-panel
Affordable performance, compact design. 10x oversampling, USB connectivity, 16 automated measurements.
Front-panel USB host port: Quickly store and transfer your waveforms and settings.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com