

**Product Name :**  
Cryo Vial Internal Thread

**Product Code :**  
NLAB-EDUSCHOOLEQ128001



**Description :**

Cryo Vial Internal Thread

**Technical Specification :**

Adding to the existing range of cryo vials is this vial which has internal thread in it.

The cap fitted with an O-ring screw into the vial and makes it a leak proof Cryo Vial.

These vials are provided with white, prominent graduations and a long marking area for sample identification.

Cryo Vial Internal Thread Capacity 2.0ml

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003