

**Product Name :**  
Cryo Vial

**Product Code :**  
NLAB-EDUSCHOOLEQ128002



**Description :**

Cryo Vial

**Technical Specification :**

Cryo vials are normally used for storage of biological material, human & animal cells.

These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification.

The cap of these vials is provided with a seal ring inside to make it completely leakproof.

The longitudinal grooves in the base give a non twisting feature to these vials when placed in a cryo rack.

Cryo Vial Capacity 1.0ml.

Cryo Vial Capacity 1.8ml.

Cryo Vial Capacity 4.5ml.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003