

Product Name :
Cryostat Microtome

Product Code :
NLAB-EDUSCHOOLEQ29011



Description :

Cryostat Microtome

Technical Specification :

Most versatile section cutting machine for quick sectioning of tissues.

Instantly freezes on built-in rapid freeze stage for biological, veterinary and medical diagnostic biopsies etc.

The Microtome is completely made of brass & stainless steel and is placed in stainless steel cold chamber with an observation glass window.

Provided with a hermetically sealed compressor unit with 'A' grade components to achieve quick cooling of up to -500C.

It is provided with feed range from 0-40 microns in steps of 2 microns each.

Accessories supplied with equipment.

Razor-120mm with back & handle-one

Object holder-set of three

Honing Stone-one

Oilcan with low temperature-one

Dust cover-one

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003