

Product Name : Curing Light	Product Code : NLABMEDICAL-M32100019
Description : Curing Light India, Curing Light manufacturer India, suppliers, Medical Curing Light India, Hospital Curing Light India	
Technical Specification : Power source: 230V~115V or 115V~(on request), 50/60Hz Weight: 2.4 Kgs. Dimension: W200 x D140 x H80(mm) Curing Light Function: Complete curing restoration with a 75 watt, 12V, halogen lamp Light guide with 360 degree rotation and can be autoclavable 3 seconds interval to turn on or off the lamp controlled by computer to protect and lengthen use life of lamp Individual curing time selection with beep tone signal Automatically turning off fan after cool down period Provides fuse and thermal cut-off switch to prevent overheating	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003