

Product Name :
Cut Section Model PF Sliding Mesh Gear Box

Product Code :
NLAB-VOCATIONALTVET26005



Description :

Cut Section Model PF Sliding Mesh Gear Box

Technical Specification :

Orientation

Material Mild Steel

Speed 750-1500rpm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003