

<b>Product Name :</b> Cylinder Trolley	<b>Product Code :</b> NLABMEDICAL-M2730014
<b>Description :</b>	
Cylinder Trolley	
<b>Technical Specification :</b>	
Color: White	
Cylinder Mounted: 150 mm Castors	
Made up of CRCA pipes mounted.	
Pretreated and epoxy powder coated.	
Tubular ms frame work fitted with two wheel, 10 cm dia.	
Finish in epoxy powder coated.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003