

Product Name :
Daniel Cell

Product Code :
NLAB-EDUSCHOOLEQ116002



Description :

Daniel Cell

Technical Specification :

Daniel Cell, Copper Vessel.

The zinc rod and copper vessel are each provided with a 4mm

socket terminal.

125 x 75mm dia.

Heavy-gauge copper outer vessel with an internal perforated shelf for

copper sulphate crystals, a porous pot and amalgamated zinc rod fitted with supporting cap

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003