


**Product Name :**  
Datalogger

**Product Code :**  
NLAB-VOCATIONALTVET28001



**Description :**

Datalogger

**Technical Specification :**

Sensors included: voltage, current, temperature, luminosity, magnetic field, absolute pressure of gases

Graphic Liquid-Crystal Display (LCD): 128 x 64 pixels

5-keys Joystick

Number of analog inputs for internal sensors: 4

Number of digital inputs: 3

RAM: 2 MB

Number of analog outputs: 2

Internal signal generator

Max. sampling frequency per channel: 1 MHz

It can be powered via USB when connected to a PC

USB Full Speed interface (12 Mbps) 2.0

It can be powered by an external USB power supply

Conversion resolution: 12 bit

Dimensions: 158 x 93 x 30 mm

Number of analog inputs for external sensors: 4

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003