

Product Name :
D.C Ammeter

Product Code :
ES741



Description :

D.C Ammeter

Technical Specification :

It is a laboratory apparatus for measurement of dc current in electric circuit for use in physics lab experiments
It is a analog desk type meter for the range 0-10 amperes with minimum error
Range: 0-10 A.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003