

**Product Name :**  
D.C Ammeters

**Product Code :**  
ES345



**Description :**

D.C Ammeters

**Technical Specification :**

It is a laboratory apparatus for measurement of dc current in electric circuit for use in physics lab experiments  
It is a analog desk type meter for the range 0-5 amperes with minimum error  
Range: 0-5 A.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003