

**Product Name :**  
D.C Voltmeter

**Product Code :**  
ES766



**Description :**

D.C Voltmeter

**Technical Specification :**

Our meters are mounted in black plastic cases.  
Meters have a scale length of 75mm, offer +/- 2.5% accuracy and come with 5-way binding post and external calibrating screw.  
Meters are not protected against prolonged or significant overload.  
Stackable for storage.  
Range: 0-300V.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003