

**Product Name :**  
Decade Inductance Box

**Product Code :**  
NLAB-TESTINGLABM17002



**Description :**

Decade Inductance Box

**Technical Specification :**

- Range : 10 $\mu$ H to 111,1mH ( 10 $\mu$ H per step )
- Accuracy : 5% inductors used throughout
- With four decades of inductance
- Max. rating current : 10mA DC or AC Current
- Internal residual inductance : Approx. 0.5 $\mu$ H
- Accurate and reliable performance
- Slide switches to allow the user to simply add or subtract for desired value
- Terminals with multi way binding posts, one to switch shield case
- ABS plastic housing case, rugged components

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003