

Product Name :
Decade Resistance Box 6 Dials

Product Code :
NLAB-TESTINGLABM17004



Description :

Decade Resistance Box 6 Dials

Technical Specification :

Specifications:

Metal Film Resistor Accuracy : $\pm 1\%$

Power Rating : 1Watt

Maximum Working Voltage : 500 Volts.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003