

Product Name : Deep Fat Fryer	Product Code : NLAB-VOCATIONALTVET27018
Description :	
Deep Fat Fryer	
Technical Specification :	
Automation Grade: Semi-Automatic	
Frequency (Hz): 50 Hz	
Material: Stainless Steel	
Power: 2000W	
Thermostatic Temperature Control: 130 Degree C to 190 Degree C	
Voltage: 220-240V	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003