

Product Name :
Deep Freezer

Product Code :
NLAB-SCIENTIFICEQ35003



Description :

Deep Freezer

Technical Specification :

- Microprocessor based control system.(HMI+PLC optional)
- Defrost cycle can be initiated to shut down the compressor and keep the cabinet defrost free.
- Ergonomic and easy to use.
- Fully programmable PLC and USB ports for PC connectivity.
- RS 232/485 connectivity
- CFC/HCFC-free refrigerants maximize cooling and minimize environmental impact
- chart recorder for better temperature monitoring (optional)
- Corrosion resistant interior and exterior
- Incident monitor provides visual and audible warnings of system abnormalities for increased sample

protection.

- Digital display and direct settings for maintaining desired operations.
- Castors for ease of movement
- Lockable hinged doors

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003