

**Product Name :**  
Deep Freezer

**Product Code :**  
NLAB-SCIENTIFICEQ35003



**Description :**

Deep Freezer

**Technical Specification :**

- Microprocessor based control system.(HMI+PLC optional)
- Defrost cycle can be initiated to shut down the compressor and keep the cabinet defrost free.
- Ergonomic and easy to use.
- Fully programmable PLC and USB ports for PC connectivity.
- RS 232/485 connectivity
- CFC/HCFC-free refrigerants maximize cooling and minimize environmental impact
- chart recorder for better temperature monitoring (optional)
- Corrosion resistant interior and exterior
- Incident monitor provides visual and audible warnings of system abnormalities for increased sample

protection.

- Digital display and direct settings for maintaining desired operations.
- Castors for ease of movement
- Lockable hinged doors

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003