

Product Name :
Deep Throat Swaging & Cosing Machine

Product Code :
NLAB-MACHINELABWK23001



Description :

Deep Throat Swaging & Cosing Machine

Technical Specification :

Deep Throat Swaging & Cosing Machine .

Maximum Capacity in M. S.

~~18~~ For hand operated SWG.

~~11~~ For motorised SWG.

~~16~~ Maximum depth of gap to gauge
inch

~~20~~ Working centres of spindles inch

~~2~~ H.P. required

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003