

Product Name :
Deflection Magnetometer

Product Code :
NLAB-SCIENTIFICEQ11004



Description :

Deflection Magnetometer

Technical Specification :

- Comprising a stout bakelite box, 100mm dia, with aluminum dial graduated in degrees from 0 to 360, anti parallax mirror slot & replaceable steel pivot.
- The short magnetic needle has a jewel bearing & a light aluminum pointer designed for precise reading & maximum damping.
- One meter long single piece wooden nicely polished base with embedded two half meter scale reading from the center of compass box i.e. zero at center of both the scales.
- Complete with compass box & leveling screw.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003