

Product Name : Demonstration Lysimeter	Product Code : NLAB-TECHNICALAB27004
Description :	
Demonstration Lysimeter	
Technical Specification :	
<p>Base and inner disc to support soil filled container and plant, resting on a hydraulic sensing device located in the base that is connected to a graduated water column mounted above the lysimeter; Complete with spare plant containers, calibration weights,user instructions and manual.</p> <p>The Demonstration Lysimeter consists of three 300mm diameter containers in which may be placed any soil type and a variety of crop types grown. Each container can then in turn be placed on a hydraulically mounted plate which is used to monitor system weight changes arising from evapotranspiration. Despite its simplicity, the apparatus has a high degree of sensitivity and may be used for student project work as well as for demonstration purposes.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003