


**Product Name :**  
Demonstration of Standing Waves on a Rope

**Product Code :**  
NLAB-SCHOOLABS15012



**Description :**

Demonstration of Standing Waves on a Rope

**Technical Specification :**

To study standing wave of transverse waves in rope. Approximately 120 cm long. AC generator box 220 VAC 50 Hz. A rope looped through a pulley at one end is attached to the frequency generator box. The other side is tied to a spring scale. Can try to find the speed of the waves in the rope and the frequency from the wave source. Well packed in wooden box

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003