


Product Name :
Demonstration of Standing Waves on a Rope

Product Code :
NLAB-SCHOOLABS15012



Description :

Demonstration of Standing Waves on a Rope

Technical Specification :

To study standing wave of transverse waves in rope. Approximately 120 cm long. AC generator box 220 VAC 50 Hz. A rope looped through a pulley at one end is attached to the frequency generator box. The other side is tied to a spring scale. Can try to find the speed of the waves in the rope and the frequency from the wave source. Well packed in wooden box

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003