

**Product Name :**  
Demonstration Set for Flowmeters

**Product Code :**  
NLAB-EDUCATIONTVET14007



**Description :**

Demonstration Set for Flowmeters

**Technical Specification :**

Allows:

- \* Directly compare the measurements of Venturi flowmeters, variable section and orifice
- \* Calibrate each flow meter using the volumetric measuring tank incorporated in the hydraulic bench
- \* Compare pressure drops across each flow meter

The flowmeters are connected in series to directly compare the measured values. The flow variation is controlled by a valve

Pressure measurement points are integrated in the flowmeters so that the pressure drops characteristic of each flowmeter can be measured.

These are measured by 8 manometric tubes with air purge valve and a hand pump allows for back pressure

Technical characteristics:

Pressure gauges: 0 to 440 mm

Number of pressure tubes 8

Diameter of the hole in the plate: 20 mm

Variable area flowmeter 2 to 20 liters / min

Venturi dimensions:

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003