

Product Name :
Demonstration Set Of Bernoulli's Theorem

Product Code :
NLAB-EDUCATIONTVET14009



Description :

Demonstration Set Of Bernoulli's Theorem

Technical Specification :

The set consists of a Venturi tube made of transparent acrylic material. The pressure is measured at 8 places along the tube with an 8-tube manometric panel.

The flow rate is regulated with a valve. Manometric tube length 300 mm, the tubes can be pressurized with a hand pump. Venturi tube diameters: 25 mm at the inlet and 10 mm at the throat

14 ° upstream cone, 21 ° downstream cone

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003